



◎ VesTraveler® Design Features

- **Semi-Compliant Balloon Catheter:** designed specifically for Pre-dilatation and Drug/Stent delivery
- **Excellent Crossability:** Unique taper-tip design, easy using in small vessel.
- **Great Trace/Travelability:** Thin wall balloon with clear marker band-lowest profile even coated with drug or crimped with stents, easy traveling and tracing inside vessel.
- **Precisely controlled balloon compliance and length expansion:** drugs and stents can be delivered as intended.
- **Maximum Pushability/ Deliverability:** TranSmooth™ Technology provides maximum pushability in most torque vessels.




◎ VesTraveler® Compliance Chart



Dongguang TT Medical, Inc.

Guangdong, PR China | Info@ttmedicalinc.com
 +86 0769-22215363 | www.ttmedicalinc.com



Technical Information

VesTraveler® Semi-Compliant PTCA Balloon Catheter	
Platform	Rapid Exchange(Rx)
Balloon Material	PA
Marker Material	Platinum-Iridium
Tip Length	3 ± 1 mm
Distal Shaft Diameter	0.89 ± 0.05 mm
Proximal Shaft Diameter	0.67 ± 0.05 mm
Compliance	Semi-compliant
Nominal Pressure(NP)	8 atm
Rated Burst Pressure(RBP)	16 atm
Catheter Shaft Length	140 ± 5 cm
Maximum Guide Wire Compatibility	0.014 inch
Shelf Life	24 months
Tip Entry Profile	0.018 inch
Balloon Crossing Profile	2.00-2.25 : 0.80 mm 2.50-2.75 : 0.85 mm 3.00-3.50 : 0.90 mm 3.75 : 0.95 mm 4.00 : 1.00 mm
Minimum Guiding Catheter ID	6F / 5F

Order Information

VesTraveler® Semi-Compliant PTCA Balloon Catheter Order Information								
Balloon Diameter (mm)	Balloon Length (mm)							
	10	14	18	22	26	30	34	38
2.00	SD2.00X10	SD2.00X14	SD2.00X18	SD2.00X22	SD2.00X26	SD2.00X30	SD2.00X34	SD2.00X38
2.25	SD2.25X10	SD2.25X14	SD2.25X18	SD2.25X22	SD2.25X26	SD2.25X30	SD2.25X34	SD2.25X38
2.50	SD2.50X10	SD2.50X14	SD2.50X18	SD2.50X22	SD2.50X26	SD2.50X30	SD2.50X34	SD2.50X38
2.75	SD2.75X10	SD2.75X14	SD2.75X18	SD2.75X22	SD2.75X26	SD2.75X30	SD2.75X34	SD2.75X38
3.00	SD3.00X10	SD3.00X14	SD3.00X18	SD3.00X22	SD3.00X26	SD3.00X30	SD3.00X34	SD3.00X38
3.25	SD3.25X10	SD3.25X14	SD3.25X18	SD3.25X22	SD3.25X26	SD3.25X30	SD3.25X34	SD3.25X38
3.50	SD3.50X10	SD3.50X14	SD3.50X18	SD3.50X22	SD3.50X26	SD3.50X30	SD3.50X34	SD3.50X38
3.75	SD3.75X10	SD3.75X14	SD3.75X18	SD3.75X22	SD3.75X26	SD3.75X30	SD3.75X34	SD3.75X38
4.00	SD4.00X10	SD4.00X14	SD4.00X18	SD4.00X22	SD4.00X26	SD4.00X30	SD4.00X34	SD4.00X38

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

VesTraveler® is a trademark or registered trademark of TT Medical. All cited trademarks are the property of their respective owners.

Dongguang TT Medical, Inc.

Guangdong, PR China Info@ttmedicalinc.com

+86 0769-22215363 www.ttmedicalinc.com

